Duroval[®] K_B 8,2

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
Reagent to determine the base capacity in water

The Duroval $K_{B\ 8.2}$ method offers a simple and rapid way to determine base capacity in water of up to pH 8.2 without additional equipment.

Each kit contains:

- 1 graduated measuring tube
- 1 stopper to close the measuring tube
- 1 dropper bottle indicator
- 1 bottle containing titration solution
- 1 test pipette, graduated in mmol/l.

The indicator dissolves in water with a pH value less than 8.2 to give a colourless solution. By adding the titration solution, the base capacity is neutralized, and the colour turns from colourless to red-violet. The amount of titration solution added to cause the colour to change from colourless to red-violet using the test pipette indicates the base capacity of the water sample. The base capacity value can be read off on the graduated test pipette.

Directions:

1. Rinse out the test tube with the water sample and fill up to the 5 ml mark (fig.1).
2. Add 1 drop of indicator and dissolve by swirling (fig. 2). If the sample turns red-violet, no base capacity is present.
3. If the sample remains colourless, close the test tube with the stopper (fig. 3).
4. Aspirate titration solution from the storage bottle with the graduated pipette (fig. 4 a) until the lower rim of the black piston ring is at the same level as the 0-mark (uppermost mark on the graduated scale) (fig. 4 b). Fill the graduated pipette with titration solution up to the piston ring (without air bubbles) to guarantee an accurate titration result. Air bubbles can be easily eliminated by depressing the plunger fully with its contents into the storage bottle, and then by refilling to the 0-mark. The tip of the pipette must be kept in the titration solution even when injecting it back.
5. With a slight twist firmly insert the pipette tip into the orifice of the test tube stopper (fig. 5 a). Whilst shaking (fig. 5 b), depress the pipette plunger very slowly (fig. 5 c) until the colour of the water sample turns from colourless to red-violet (fig. 6 a).
6. Read the degree of base capacity taken at the lower rim of the black piston ring on the graduated pipette (fig. 6 b) and take note. If the contents of one pipette are insufficient to effect the colour change, refill the pipette according to 4 with fresh titration solution up to the 0-mark as described and continue the titration in the same water sample according to 5 and 6.

NOTE!

Close the container tightly after drawing the titration solution.

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