# Duroval<sup>®</sup> K<sub>S</sub> 4,3

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
## Reagent to determine the acid capacity in water

The Duroval  $K_{S\ 4,3}$  method offers a simple and rapid method to determine acid capacity in water.

Each kit contains:

- 1 graduated measuring tube
- 1 stopper to close the measuring tube
- 1 dropper bottle indicator
- 1 bottle containing titration solution
- 1 test pipette, graduated in mmol/l.

The indicator disperses in water with a pH value greater than 4.3 to blue-green. By adding the titration solution, the acid capacity is neutralized, and the colour turns from blue-green to orange. The amount of titration solution added to cause the colour to change from blue-green to orange using the test pipette indicates the acid capacity of the water sample. The acid capacity value can be read off on the graduated test pipette.

### Directions:

1. Rinse out the test tube with the water sample and fill up to the 5 ml mark (fig.1).
2. Add 1 drop of indicator and disperse by swirling (fig. 2). If the sample turns to orange, no acid capacity is present.
3. If the sample turns to blue-green, close the test tube with the stopper (fig. 3).
4. Aspirate titration solution from the storage bottle with the graduated pipette (fig. 4 a) until the lower rim of the black piston ring is at the same level as the 0-mark (uppermost mark on the graduated scale) (fig. 4 b). Fill the graduated pipette with titration solution up to the piston ring (without air bubbles) to guarantee an accurate titration result. Air bubbles can be easily eliminated by depressing the plunger fully with its contents into the storage bottle, and then by refilling to the 0-mark. The tip of the pipette must be kept in the titration solution even when injecting it back.
5. With a slight twist firmly insert the pipette tip into the orifice of the test tube stopper (fig. 5 a). Whilst shaking (fig. 5 b), depress the pipette plunger very slowly (fig. 5 c) until the colour of the water sample turns from blue-green to orange (fig. 6 a).
6. Read the degree of acid capacity taken at the lower rim of the black piston ring on the graduated pipette (fig. 6 b) and take note. If the contents of one pipette are insufficient to effect the colour change, refill the pipette according to 4 with fresh titration solution up to the 0-mark as described and continue the titration in the same water sample according to 5 and 6.

Note! When determining the acid capacity, allow for the quantity of titration solution previously used.

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