testoval[®]

ALUMINIUM - TEST

Al³⁺ gelöst

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Directions for use

Remove the test cell from the comparator and fill it up to the mark with water sample. Add 1 drop of reagent A and mix by stirring with the red spoon. After 1 minute add 10 drops of reagent B to the water sample and mix. Reinsert the test cell into the comparator and after an interval of 5 minutes hold the comparator against a bright background or light source and match the colour intensity of the central field (test solution) with one of the colours of the peripheral colour fields. Read the value in ppm Al, which is printed under the matching colour field.

Should the colour in the central field be more intensive than any of the colours in the peripheral fields, the water sample should be diluted 1:10. For this purpose fill the graduated plastic tube up to the 1 ml-mark with sample water and dilute with aluminium-free water (e.g. distilled water) to the 10 ml-mark. Transfer the diluted sample into the test cell, which has been rinsed with aluminium-free water and proceed as described above. The value obtained must in this case be multiplied by 10.

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