


# Safety data sheet according to 1907/2006/EC

Printing date 06.08.2021

Version number 2 (replaces version 1)

Revision: 06.08.2021



- Signal word Void
- · Hazard statements Void
- · 2.3 Other hazards
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- vPvB: Not applicable.

#### SECTION 3: Composition/information on ingredients

- · 3.2 Mixtures
- · Description:

Mixture of substances listed below with nonhazardous additions according to Regulation (EC) No 1272/2008. Sodium acetate CAS 6131-90-4

Water CAS 7732-18-5

- · Dangerous components: Void
- · SVHC Not applicable.
- Additional information: For the wording of the listed hazard phrases refer to section 16.

#### SECTION 4: First aid measures

• 4.1 Description of first aid measures

- · General information: Immediately remove any clothing soiled by the product.
- After inhalation: Supply fresh air; consult doctor in case of complaints.

(Contd. on page 2)

DE\_EN



## Printing date 06.08.2021 Ve

according to 1907/2006/EC Version number 2 (replaces version 1)

Safety data sheet

Revision: 06.08.2021

## Trade name: Testomat Chlor Reagenz AF

	(Contd. of page 1
After skin contact:	ίται το σ
Immediately wash with water and soap and rinse thoroughly.	
If skin irritation continues, consult a doctor.	
After eye contact:	
Rinse opened eye for several minutes under running water. It <b>After swallowing:</b>	f symptoms persist, consult a doctor.
A person vomiting while laying on their back should be turned	d onto their side
Rinse out mouth and then drink plenty of water.	
Seek medical treatment.	
4.2 Most important symptoms and effects, both acute an	d delaved
No further relevant information available.	
4.3 Indication of any immediate medical attention and sp	ecial treatment needed
No further relevant information available.	
SECTION 5: Firefighting measures	
5.1 Extinguishing media	
Suitable extinguishing agents:	
CO2, powder or water spray. Fight larger fires with water spray	
5.2 Special hazards arising from the substance or mixtur	re la
In case of fire, the following can be released:	
Acetic acid	
5.3 Advice for firefighters	
<b>Protective equipment:</b> Wear self-contained respiratory protection	ective device.

### SECTION 6: Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures Wear personal protection equipment.
- 6.2 Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- 6.3 Methods and material for containment and cleaning up:

Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust). Dispose of the material collected according to regulations. Clean the affected area carefully; suitable cleaners are: Warm water

**6.4 Reference to other sections** See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment. See Section 13 for disposal information.

### SECTION 7: Handling and storage

- · 7.1 Precautions for safe handling No special measures required.
- · Information about fire and explosion protection: No special measures required.
- · 7.2 Conditions for safe storage, including any incompatibilities
- · Storage:
- · Requirements to be met by storerooms and receptacles: Store only in the original receptacle.
- · Information about storage in one common storage facility: Store away from foodstuffs.
- · Further information about storage conditions:
- Keep container tightly sealed.
- Protect from heat and direct sunlight.

Store receptacle in a well ventilated area.

· Recommended storage temperature: 15 - 25 °C

(Contd. on page 3)



## Safety data sheet according to 1907/2006/EC

Printing date 06.08.2021

Version number 2 (replaces version 1)

Revision: 06.08.2021

(Contd. of page 2)

### Trade name: Testomat Chlor Reagenz AF

#### • Storage class:

- Storage class 12: Non-combustible liquids that cannot be assigned to any other storage class (TRGS 510)
- 7.3 Specific end use(s) No further relevant information available.
- SECTION 8: Exposure controls/personal protection · 8.1 Control parameters · Ingredients with limit values that require monitoring at the workplace: The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace. Additional information: The lists valid during the making were used as basis. · 8.2 Exposure controls · Appropriate engineering controls No further data; see item 7. · Individual protection measures, such as personal protective equipment General protective and hygienic measures: The usual precautionary measures are to be adhered to when handling chemicals. Keep away from foodstuffs, beverages and feed. Immediately remove all soiled and contaminated clothing Wash hands before breaks and at the end of work. Do not eat, drink, smoke or sniff while working. Respiratory protection: Use suitable respiratory protective device when aerosol or mist is formed. Filter: Type P2 Hand protection Wear gloves according to EN 374. The glove material has to be impermeable and resistant to the product/ the substance/ the preparation. Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation Check protective gloves prior to each use for their proper condition. Preventive skin protection by use of skin-protecting agents is recommended. After use of gloves apply skin-cleaning agents and skin cosmetics. Material of gloves The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.
  - Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

• For the permanent contact gloves made of the following materials are suitable: Nitrile rubber, NBR

Recommended thickness of the material:  $\geq 0.12$  mm Value for the permeation: Level = 6 (> 480 min)

• As protection from splashes gloves made of the following materials are suitable: Nitrile rubber, NBR

Recommended thickness of the material:  $\geq 0.12$  mm

Value for the permeation: Level = 6 (> 480 min)

- Eye/face protection Goggles according to EN 166 recommended during refilling
- Body protection: Protective work clothing

#### SECTION 9: Physical and chemical properties

· 9.1 Information on basic physical and chemical properties

<sup>.</sup> General Information

Physical state
Colour:

Fluid Colourless

(Contd. on page 4)

DE\_EN



Printing date 06.08.2021

## Safety data sheet according to 1907/2006/EC

Version number 2 (replaces version 1)

Revision: 06.08.2021

## Trade name: Testomat Chlor Reagenz AF

	(Contd. of pag
Odour:	Characteristic
Odour threshold:	Not determined.
Melting point/freezing point:	Undetermined.
Boiling point or initial boiling point and	
boiling range	Undetermined.
Flammability	Not applicable.
Lower and upper explosion limit	
Lower:	Not determined.
Upper:	Not determined.
Flash point:	Undetermined.
Auto-ignition temperature:	Product is not selfigniting.
Decomposition temperature:	Not determined.
pH (100 g/l) at 20 °C	5.8
Viscosity:	0.0
	Not determined.
Kinematic viscosity	
Dynamic:	Not determined.
Solubility	Fully, we'r eilele
water:	Fully miscible.
Partition coefficient n-octanol/water (log	
value)	Not determined.
Vapour pressure:	Not determined.
Density and/or relative density	
Density at 20 °C:	1.06 g/cm³
Relative density	Not determined.
Vapour density	Not determined.
and environment, and on safety. Ignition temperature:	Not determined.
Explosive properties:	
	Product does not present an explosion hazard.
Change in condition	
	Product does not present an explosion hazard. Not determined.
Change in condition Evaporation rate Information with regard to physical haz	Not determined.
Change in condition Evaporation rate Information with regard to physical haz classes	Not determined. ard
Change in condition Evaporation rate Information with regard to physical haz classes Explosives	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives	Not determined. ard
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols	Not determined. ard Void Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases	Not determined. ard Void Void Void Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure	Not determined. ard Void Void Void Void Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids	Not determined. ard Void Void Void Void Void Void Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids	Not determined. ard Void Void Void Void Void Void Void Voi
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures	Not determined. ard Void Void Void Void Void Void Void Void Void Void Void Void Void Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids Pyrophoric solids	Not determined. ard Void Vo
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids Pyrophoric solids Self-heating substances and mixtures	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids Pyrophoric solids Self-heating substances and mixtures Substances and mixtures, which emit	Not determined. ard Void Vo
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids Pyrophoric solids Self-heating substances and mixtures Substances and mixtures, which emit flammable gases in contact with water	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids Pyrophoric solids Self-heating substances and mixtures Substances and mixtures, which emit flammable gases in contact with water Oxidising liquids	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids Pyrophoric solids Self-heating substances and mixtures Substances and mixtures, which emit flammable gases in contact with water Oxidising liquids Oxidising solids	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz classes Explosives Flammable gases Aerosols Oxidising gases Gases under pressure Flammable liquids Flammable solids Self-reactive substances and mixtures Pyrophoric liquids Pyrophoric solids Self-heating substances and mixtures Substances and mixtures, which emit flammable gases in contact with water Oxidising liquids Oxidising solids Organic peroxides	Not determined. ard Void
Change in condition Evaporation rate Information with regard to physical haz	Not determined. ard Void

(Contd. on page 5)



#### Printing date 06.08.2021

## Safety data sheet according to 1907/2006/EC

Version number 2 (replaces version 1)

Revision: 06.08.2021

### Trade name: Testomat Chlor Reagenz AF

(Contd. of page 4)

#### SECTION 10: Stability and reactivity

- · 10.1 Reactivity No further relevant information available.
- 10.2 Chemical stability
- · Thermal decomposition / conditions to be avoided:
- No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid** No further relevant information available.
- 10.5 Incompatible materials: No further relevant information available.
- 10.6 Hazardous decomposition products: Acetic acid

#### **SECTION 11: Toxicological information**

- 11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008
- · Acute toxicity Based on available data, the classification criteria are not met.
- Skin corrosion/irritation Based on available data, the classification criteria are not met.
- · Serious eye damage/irritation Based on available data, the classification criteria are not met.
- · Respiratory or skin sensitisation Based on available data, the classification criteria are not met.
- Germ cell mutagenicity Based on available data, the classification criteria are not met.
- · Carcinogenicity Based on available data, the classification criteria are not met.
- Reproductive toxicity Based on available data, the classification criteria are not met.
- STOT-single exposure Based on available data, the classification criteria are not met.
- STOT-repeated exposure Based on available data, the classification criteria are not met.
- Aspiration hazard Based on available data, the classification criteria are not met.
- 11.2 Information on other hazards
- · Endocrine disrupting properties

None of the ingredients is listed.

#### **SECTION 12: Ecological information**

#### · 12.1 Toxicity

- · Aquatic toxicity: No further relevant information available.
- 12.2 Persistence and degradability No further relevant information available.
- · 12.3 Bioaccumulative potential No further relevant information available.
- · 12.4 Mobility in soil No further relevant information available.
- 12.5 Results of PBT and vPvB assessment
- **PBT:** Not applicable.
- · vPvB: Not applicable.
- 12.6 Endocrine disrupting properties
- The product does not contain substances with endocrine disrupting properties.
- 12.7 Other adverse effects
- · Additional ecological information:
- · General notes:

Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

#### **SECTION 13: Disposal considerations**

· 13.1 Waste treatment methods

· Recommendation Disposal must be made according to official regulations.

(Contd. on page 6)



## Safety data sheet according to 1907/2006/EC

Printing date 06.08.2021

Version number 2 (replaces version 1)

Revision: 06.08.2021

### Trade name: Testomat Chlor Reagenz AF

(Contd. of page 5)

- Uncleaned packaging:
- Recommendation:
- Packagings that may not be cleansed are to be disposed of in the same manner as the product. Disposal must be made according to official regulations.
- · Recommended cleansing agents: Water, if necessary together with cleansing agents.

#### SECTION 14: Transport information · 14.1 UN number or ID number · ADR. ADN. IMDG. IATA Void 14.2 UN proper shipping name · ADR, ADN, IMDG, IATA Void · 14.3 Transport hazard class(es) · ADR, ADN, IMDG, IATA · Class Void 14.4 Packing group · ADR, IMDG, IATA Void · 14.5 Environmental hazards: · Marine pollutant: No · 14.6 Special precautions for user Not applicable. · 14.7 Maritime transport in bulk according to IMO instruments Not applicable. Not dangerous according to the above Transport/Additional information: specifications. **UN "Model Regulation":** Void

## SECTION 15: Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture
- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I None of the ingredients is listed.
- DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II

None of the ingredients is listed.

· REGULATION (EU) 2019/1148

 Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))


None of the ingredients is listed.

#### · Annex II - REPORTABLE EXPLOSIVES PRECURSORS

None of the ingredients is listed.

#### National regulations:

#### · Technical instructions (air):



(Contd. on page 7)

DE\_EN-



Printing date 06.08.2021

## Safety data sheet according to 1907/2006/EC

Version number 2 (replaces version 1)

Revision: 06.08.2021

### Trade name: Testomat Chlor Reagenz AF

(Contd. of page 6)

• Waterhazard class: Water hazard class 1 (Self-assessment): slightly hazardous for water. • 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

#### **SECTION 16: Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· **Department issuing SDS:** product safety department

Date of previous version: 19.05.2020

- · Version number of previous version: 1
- · Abbreviations and acronvms: RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail) ICAO: International Civil Aviation Organisation ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road) IMDG: International Maritime Code for Dangerous Goods IATA: International Air Transport Association GHS: Globally Harmonised System of Classification and Labelling of Chemicals EINECS: European Inventory of Existing Commercial Chemical Substances ELINCS: European List of Notified Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society) PBT: Persistent, Bioaccumulative and Toxic SVHC: Substances of Very High Concern vPvB: very Persistent and very Bioaccumulative \*\* Data compared to the previous version altered. DE EN