

Version number 10 (replaces version 9)

Revision: 07.12.2021

SECTION 1: Identification of the substance/mixture undertaking	e and of the company/
· 1.1 Product identifier	
<ul> <li>Trade name: <u>Testoval Eisen Test 0-10 mg/l B</u></li> <li>CAS Number: 50-81-7</li> <li>EC number: 200-066-2</li> <li>UFI: 0D1S-ETMY-UC0F-5CQP</li> </ul>	
<ul> <li>• 1.2 Relevant identified uses of the substance or mixture and uses No further relevant information available.</li> <li>• Application of the substance / the preparation: Reagent for analysis EuPCS: PC-TEC-19 Reagents and laboratory chemical</li> </ul>	s advised against
<ul> <li>1.3 Details of the supplier of the safety data sheet</li> <li>Manufacturer/Supplier: Gebrüder Heyl Analysentechnik GmbH &amp; Co. KG Orleansstraße 75 b D-31135 Hildesheim</li> <li>Further information obtainable from: product safety department</li> </ul>	Phone +49 (0) 5121 2893390 Fax +49 (0) 5121 2893367 E-mail info@heylanalysis.de Internet www.heylanalysis.de
• <b>1.4 Emergency telephone number:</b> Giftinformationszentrum Nord Phone +49 (0) 551 19240	
SECTION 2: Hazards identification	
<ul> <li>2.1 Classification of the substance or mixture</li> <li>Classification according to Regulation (EC) No 1272/2008</li> <li>The substance is not classified, according to the CLP regulation.</li> </ul>	
<ul> <li>2.2 Label elements</li> <li>Labelling according to Regulation (EC) No 1272/2008 Void</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>Labelling of packages where the contents do not exceed 125 ml</li> <li>Hazard pictograms Void</li> <li>Signal word Void</li> <li>Signal word Void</li> <li>Signal word Void</li> <li>Hazard statements Void</li> <li>Signal word Void</li> <li>Results of PBT and vPvB assessment</li> <li>PBT: Not applicable.</li> <li>vPvB: Not applicable.</li> </ul>	
SECTION 3: Composition/information on ingredients	
<ul> <li>• 3.1 Substances</li> <li>• CAS No. Description CAS: 50-81-7 ascorbic acid</li> <li>• Identification number(s)</li> </ul>	

• SVHC Not applicable.

\_\_\_\_\_ DE\_EN (Contd. on page 2)



# Safety data sheet according to 1907/2006/EC

Version number 10 (replaces version 9)

Revision: 07.12.2021

#### Trade name: Testoval Eisen Test 0-10 mg/l B

(Contd. of page 1)

## **SECTION 4: First aid measures**

- · 4.1 Description of first aid measures
- General information: Immediately remove any clothing soiled by the product.
- After inhalation: Supply fresh air; consult doctor in case of complaints.
- After skin contact:
- Immediately rinse with water.

If skin irritation continues, consult a doctor.

• After eye contact:

*Rinse opened eye for several minutes under running water. If symptoms persist, consult a doctor. After swallowing:* 

A person vomiting while laying on their back should be turned onto their side. Rinse out mouth and then drink plenty of water.

- Seek medical treatment.
- **4.2 Most important symptoms and effects, both acute and delayed** No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed** No further relevant information available.

## **SECTION 5: Firefighting measures**

- · 5.1 Extinguishing media
- · Suitable extinguishing agents: Use fire extinguishing methods suitable to surrounding conditions.
- · For safety reasons unsuitable extinguishing agents: Water with full jet
- 5.2 Special hazards arising from the substance or mixture
- Formation of toxic gases is possible during heating or in case of fire.
- 5.3 Advice for firefighters
- · Protective equipment: Wear self-contained respiratory protective device.

## SECTION 6: Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures Wear personal protection equipment. Ensure adequate ventilation Avoid formation of dust.
   6.2 Environmental precautions: Do not allow to enter sewers/ surface or ground water.
- 6.3 Methods and material for containment and cleaning up:
- Pick up mechanically.

Dispose contaminated material as waste according to item 13. Clean the affected area carefully; suitable cleaners are:

- Warm water
- **6.4 Reference to other sections** See Section 7 for information on safe handling. See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

## SECTION 7: Handling and storage

• **7.1 Precautions for safe handling** Ensure good ventilation/exhaustion at the workplace. Prevent formation of dust. Any unavoidable deposit of dust must be regularly removed.

· Information about fire - and explosion protection: No special measures required.

(Contd. on page 3)

DE EN-



## Safety data sheet according to 1907/2006/EC

Version number 10 (replaces version 9)

Revision: 07.12.2021

## Trade name: Testoval Eisen Test 0-10 mg/l B

(Contd. of page 2)

- · 7.2 Conditions for safe storage, including any incompatibilities
- · Storage:
- · Requirements to be met by storerooms and receptacles: Store only in the original receptacle.
- · Information about storage in one common storage facility: Store away from foodstuffs.
- Further information about storage conditions:
- Keep container tightly sealed. Protect from heat and direct sunlight. Store in dry conditions. Store receptacle in a well ventilated area.
- Recommended storage temperature: 15 25 °C
- Storage class:

Storage class 11: Combustible solids that cannot be assigned to any other storage class (TRGS 510)

• 7.3 Specific end use(s) No further relevant information available.

#### SECTION 8: Exposure controls/personal protection

- · 8.1 Control parameters
- · Ingredients with limit values that require monitoring at the workplace: Not required.
- Additional information: The lists valid during the making were used as basis.
- · 8.2 Exposure controls
- Appropriate engineering controls No further data; see item 7.
- · Individual protection measures, such as personal protective equipment
- General protective and hygienic measures:
- The usual precautionary measures are to be adhered to when handling chemicals.
- Keep away from foodstuffs, beverages and feed.
- Immediately remove all soiled and contaminated clothing

Wash hands before breaks and at the end of work.

Do not eat, drink, smoke or sniff while working.

· Respiratory protection:

Use suitable respiratory protective device when aerosol or mist is formed. Filter: Type P2

· Hand protection

Wear gloves according to EN 374.

The glove material has to be impermeable and resistant to the product/ the substance/ the preparation.

Selection of the glove material on consideration of the penetration times, rates of diffusion and the degradation

Check protective gloves prior to each use for their proper condition.

Preventive skin protection by use of skin-protecting agents is recommended.

After use of gloves apply skin-cleaning agents and skin cosmetics.

• Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer.

#### · Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

- For the permanent contact gloves made of the following materials are suitable: Nitrile rubber, NBR
- Recommended thickness of the material:  $\geq 0.12 \text{ mm}$
- Value for the permeation: Level = 6 (> 480 min)
- As protection from splashes gloves made of the following materials are suitable: Nitrile rubber, NBR

Recommended thickness of the material:  $\geq 0.12$  mm Value for the permeation: Level = 6 (> 480 min)

· Eye/face protection Not required.

(Contd. on page 4)



# Safety data sheet according to 1907/2006/EC

Version number 10 (replaces version 9)

Revision: 07.12.2021

(Contd. of page 3)

## Trade name: Testoval Eisen Test 0-10 mg/l B

· Body protection: Protective work clothing

	properties	
9.1 Information on basic physical and cher	mical properties	
General Information Physical state	Solid	
Colour:	White	
Odour:	Odourless	
Odour threshold:	Not determined.	
	191 °C	
<i>Melting point/freezing point: Boiling point or initial boiling point and</i>	191 C	
boiling range	Undetermined.	
Flammability	Product is not flammable.	
Lower and upper explosion limit		
Lower:	Not determined.	
Upper:	Not determined.	
Flash point:	Undetermined.	
Auto-ignition temperature:	Not determined.	
Decomposition temperature:	Not determined.	
pH (50 g/l) at 20 °C	ca. 2.4	
Viscosity:	<b>UU. 2.7</b>	
Kinematic viscosity	Not applicable.	
Dynamic:	Not applicable.	
Solubility		
water at 20 °C:	333 g/l	
Partition coefficient n-octanol/water (log		
value)	Not determined.	
Vapour pressure:	Not applicable.	
Density and/or relative density		
Density at 20 °C:	0.95 g/cm³	
Relative density	Not determined.	
Bulk density:	$1,700 \text{ kg/m}^3$	
Vapour density	Not applicable.	
Particle characteristics	See item 3.	
9.2 Other information	5/46	
Important information on protection of hea	11(11	
and environment, and on safety. Ignition tomporature:	380 °C	
Ignition temperature:	Not determined.	
Explosivo proportios:		
Explosive properties:	Product does not present an explosion hazard.	
Change in condition	Not applicable	
Evaporation rate	Not applicable.	
Information with regard to physical haza	ard	
classes		
Explosives	Void	
Flammable gases	Void	
Aerosols	Void	
Oxidising gases	Void	
Gases under pressure	Void	
Flammable liquids	Void	
Flammable solids	Void	
Self-reactive substances and mixtures	Void	
Pyrophoric liquids	Void	
Pyrophoric solids	Void	

page 5)



## Safety data sheet according to 1907/2006/EC

Version number 10 (replaces version 9)

Revision: 07.12.2021

## Trade name: Testoval Eisen Test 0-10 mg/l B

		(Contd. of page 4)
<ul> <li>Self-heating substances and mixtures</li> <li>Substances and mixtures, which emit</li> </ul>	Void	
flammable gases in contact with water	Void	
· Oxidising liquids	Void	
· Oxidising solids	Void	
· Organic peroxides	Void	
· Corrosive to metals	Void	
· Desensitised explosives	Void	

## SECTION 10: Stability and reactivity

· 10.1 Reactivity No further relevant information available.

- · 10.2 Chemical stability
- **Thermal decomposition / conditions to be avoided:** No decomposition if used according to specifications.
- 10.3 Possibility of hazardous reactions
- Reacts with oxidising agents.
- Reacts with light metals.
- Risk of dust explosion if enriched with fine dust in the presence of air.
- **10.4 Conditions to avoid** No further relevant information available.
- · 10.5 Incompatible materials: No further relevant information available.
- 10.6 Hazardous decomposition products: No dangerous decomposition products known.

#### **SECTION 11: Toxicological information**

11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008

· Acute toxicity Based on available data, the classification criteria are not met.

· LD/LC50 values relevant for classification:

CAS: 50-81-7 ascorbic acid

Oral LD50 11,900 mg/kg (rat)

· Skin corrosion/irritation Based on available data, the classification criteria are not met.

- Serious eye damage/irritation Based on available data, the classification criteria are not met.
- · Respiratory or skin sensitisation Based on available data, the classification criteria are not met.
- Germ cell mutagenicity Based on available data, the classification criteria are not met.

· Carcinogenicity Based on available data, the classification criteria are not met.

- Reproductive toxicity Based on available data, the classification criteria are not met.
- · STOT-single exposure Based on available data, the classification criteria are not met.
- STOT-repeated exposure Based on available data, the classification criteria are not met.
- Aspiration hazard Based on available data, the classification criteria are not met.
- 11.2 Information on other hazards

• Endocrine disrupting properties

Substance is not listed.

## SECTION 12: Ecological information

#### · 12.1 Toxicity

- Aquatic toxicity: No further relevant information available.
- 12.2 Persistence and degradability Easily biodegradable
- 12.3 Bioaccumulative potential
- Due to the distribution coefficient n-octanol/water an accumulation in organisms is not expected.
- · 12.4 Mobility in soil No further relevant information available.

(Contd. on page 6)



## Safety data sheet according to 1907/2006/EC

Version number 10 (replaces version 9)

Revision: 07.12.2021

(Contd. of page 5)

## Trade name: Testoval Eisen Test 0-10 mg/l B

· 12.5 Results of PBT and vPvB assessment

- · PBT: Not applicable.
- · vPvB: Not applicable.
- 12.6 Endocrine disrupting properties

The product does not contain substances with endocrine disrupting properties.

- 12.7 Other adverse effects
- · Additional ecological information:
- · General notes:

Water hazard class 1 (German Regulation) (Assessment by list): slightly hazardous for water Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

#### SECTION 13: Disposal considerations

· 13.1 Waste treatment methods

· Recommendation Disposal must be made according to official regulations.

- · Uncleaned packaging:
- Recommendation:

Packagings that may not be cleansed are to be disposed of in the same manner as the product. Disposal must be made according to official regulations.

SECTION 14: Transport information	tion
· 14.1 UN number or ID number · ADR, ADN, IMDG, IATA	Void
<ul> <li>14.2 UN proper shipping name</li> <li>ADR, ADN, IMDG, IATA</li> </ul>	Void
<ul> <li>14.3 Transport hazard class(es)</li> </ul>	
· ADR, ADN, IMDG, IATA · Class	Void
· 14.4 Packing group · ADR, IMDG, IATA	Void
<ul> <li>14.5 Environmental hazards:</li> <li>Marine pollutant:</li> </ul>	No
· 14.6 Special precautions for user	Not applicable.
<ul> <li>14.7 Maritime transport in bulk accordi IMO instruments</li> </ul>	ng to Not applicable.
· Transport/Additional information:	Not dangerous according to the above specifications.
· UN "Model Regulation":	Void

## SECTION 15: Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture
- · Directive 2012/18/EU
- · Named dangerous substances ANNEX I Substance is not listed.



# Safety data sheet according to 1907/2006/EC

Printing date 07.12.2021

Version number 10 (replaces version 9)

Revision: 07.12.2021

## Trade name: Testoval Eisen Test 0-10 mg/l B

(Contd. of page 6)

#### DIRECTIVE 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment – Annex II

Substance is not listed.

· REGULATION (EU) 2019/1148

 Annex I - RESTRICTED EXPLOSIVES PRECURSORS (Upper limit value for the purpose of licensing under Article 5(3))

Substance is not listed.

Annex II - REPORTABLE EXPLOSIVES PRECURSORS

Substance is not listed.

· Regulation (EC) No 273/2004 on drug precursors

Substance is not listed.

 Regulation (EC) No 111/2005 laying down rules for the monitoring of trade between the Community and third countries in drug precursors

Substance is not listed.

· National regulations:

• Waterhazard class: Water hazard class 1 (Assessment by list): slightly hazardous for water.

· 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

## **SECTION 16: Other information**

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

· Department issuing SDS: product safety department

• Date of previous version: 17.10.2019

- Version number of previous version: 9
- Abbreviations and acronyms:

RID: Règlement international concernant le transport des marchandises dangereuses par chemin de fer (Regulations Concerning the International Transport of Dangerous Goods by Rail)

ICAO: International Civil Aviation Organisation

ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)

- IMDG: International Maritime Code for Dangerous Goods
- IATA: International Air Transport Association

GHS: Globally Harmonised System of Classification and Labelling of Chemicals

EINECS: European Inventory of Existing Commercial Chemical Substances

CAS: Chemical Abstracts Service (division of the American Chemical Society)

LC50: Lethal concentration, 50 percent

LD50: Lethal dose, 50 percent

PBT: Persistent, Bioaccumulative and Toxic

SVHC: Substances of Very High Concern vPvB: very Persistent and very Bioaccumulative

\* Data compared to the previous version altered.

DE EN-